

## 兔抗 NUDT8 多克隆抗体

中文名称: 兔抗 NUDT8 多克隆抗体

英文名称: Anti-NUDT8 rabbit polyclonal antibody

别名: nudix hydrolase 8

储存: 冷冻 (-20℃)

抗原: NUDT8

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Probably mediates the hydrolysis of some nucleoside diphosphate derivatives.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NUDT8
<b>Immunogen:</b>	Synthetic peptide of human NUDT8
<b>Full name:</b>	nudix hydrolase 8
<b>SwissProt:</b>	Q8WV74
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	20-100

