

兔抗 OLA1 多克隆抗体

- 中文名称: 兔抗 OLA1 多克隆抗体
- 英文名称: Anti-OLA1 rabbit polyclonal antibody
- 别名: Obg like ATPase 1; DOC45; GBP45; GTBP9; GTPBP9; PTD004
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 抗原: OLA1
- 反应种属: Human, Mouse, Rat
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the GTPase protein family. The encoded protein interacts with breast cancer-associated gene 1 (BRCA1) and BRCA1-associated RING domain protein (BARD1), and is involved in centrosome regulation. Overexpression of this gene has been observed in multiple types of cancer and may be associated with poor survival. Pseudogenes of this gene have been defined on chromosomes 17 and 22.
Applications:	ELISA, WB, IHC
Name of antibody:	OLA1
Immunogen:	Fusion protein of human OLA1
Full name:	Obg like ATPase 1

Synonyms:	DOC45; GBP45; GTBP9; GTPBP9; PTD004
SwissProt:	Q9NTK5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human lung cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	45 kDa
WB Positive control:	A375 cell lysate
WB Recommended dilution:	500-2000



