

兔抗 TBX6 多克隆抗体

中文名称：兔抗 TBX6 多克隆抗体

英文名称：Anti-TBX6 rabbit polyclonal antibody

别名：SCDO5

抗原：TBX6

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. Knockout studies in mice indicate that this gene is important for specification of paraxial mesoderm structures.
Applications:	ELISA, IHC
Name of antibody:	TBX6
Immunogen:	Fusion protein of human TBX6
Full name:	T-box 6
Synonyms:	SCDO5

SwissProt:	O95947
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	20-100

