

兔抗 TCOF1 多克隆抗体

中文名称: 兔抗 TCOF1 多克隆抗体

英文名称: Anti-TCOF1 rabbit polyclonal antibody

别名: treacle ribosome biogenesis factor 1; TCS; MFD1; TCS1; treacle

储存: 冷冻 (-20℃)

抗原: TCOF1

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a nucleolar protein with a LIS1 homology domain. The protein is involved in ribosomal DNA gene transcription through its interaction with upstream binding factor (UBF). Mutations in this gene have been associated with Treacher Collins syndrome, a disorder which includes abnormal craniofacial development. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TCOF1

Immunogen:	Fusion protein of human TCOF1
Full name:	treacle ribosome biogenesis factor 1
Synonyms:	TCS; MFD1; TCS1; treacle
SwissProt:	Q13428
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	40-200



