

兔抗 GINS4 多克隆抗体

中文名称： 兔抗 GINS4 多克隆抗体

英文名称： Anti-GINS4 rabbit polyclonal antibody

别名： SLD5

抗原： GINS4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The yeast heterotetrameric GINS complex is made up of Sld5, Psf1 (GINS1; MIM 610608), Psf2 (GINS2; MIM 610609), and Psf3 (GINS3; MIM 610610). The formation of the GINS complex is essential for the initiation of DNA replication in yeast and <i>Xenopus</i> egg extracts (Ueno et al., 2005 [PubMed 16287864]). See GINS1 for additional information about the GINS complex.
Applications:	ELISA, IHC
Name of antibody:	GINS4
Immunogen:	Full length fusion protein
Full name:	GINS complex subunit 4 (Sld5 homolog)
Synonyms :	SLD5

SwissProt:	Q9BRT9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	25-100

