

兔抗 GJB5 多克隆抗体

中文名称: 兔抗 GJB5 多克隆抗体

英文名称: Anti-GJB5 rabbit polyclonal antibody

别名: CX31.1

储存: 冷冻 (-20℃)

抗原: GJB5

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the beta-type (group I) connexin family. The encoded protein is a gap junction protein involved in intercellular communication related to epidermal differentiation and environmental sensing. This gene has been linked to non-small cell lung cancer.
Applications:	ELISA, WB
Name of antibody:	GJB5
Immunogen:	Synthetic peptide of human GJB5
Full name:	gap junction protein, beta 5, 31.1kDa
Synonyms :	CX31.1
SwissProt:	O95377

ELISA Recommended dilution:	1000-2000
WB Predicted band size:	31 kDa
WB Positive control:	Human placenta tissue
WB Recommended dilution:	200-1000

