

兔抗 UBE4B 多克隆抗体

中文名称:	兔抗 UBE4B 多克隆抗体
英文名称:	Anti-UBE4B rabbit polyclonal antibody
别名:	E4; UFD2; HDNB1; UBOX3; UFD2A
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	UBE4B
反应种属:	Human, Mouse
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes an additional conjugation factor, E4, which is involved in multiubiquitin chain assembly. This gene is also the strongest candidate in the neuroblastoma tumor suppressor genes. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	UBE4B
Immunogen:	Synthetic peptide of human UBE4B
Full name:	ubiquitination factor E4B
Synonyms:	E4; UFD2; HDNB1; UBOX3; UFD2A
SwissProt:	O95155
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

