

## 兔抗 UBXN1 多克隆抗体

中文名称：兔抗 UBXN1 多克隆抗体

英文名称：Anti-UBXN1 rabbit polyclonal antibody

别名：UBX domain protein 1; 2B28; SAKS1; UBXD10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UBXN1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Ubiquitin-binding protein that interacts with the BRCA1-BARD1 heterodimer, and regulates its activity. Specifically binds 'Lys-6'-linked polyubiquitin chains. Interaction with autoubiquitinated BRCA1, leads to inhibit the E3 ubiquitin-protein ligase activity of the BRCA1-BARD1 heterodimer. Component of a complex required to couple deglycosylation and proteasome-mediated degradation of misfolded proteins in the endoplasmic reticulum that are retrotranslocated in the cytosol.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	UBXN1
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	UBX domain protein 1
<b>Synonyms:</b>	2B28; SAKS1; UBXD10

<b>SwissProt:</b>	Q04323
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	33 kDa
<b>WB Positive control:</b>	Mouse brain tissue , Hela and HEPG2 cell lysates
<b>WB Recommended dilution:</b>	500-2000





