

兔抗 UFC1 多克隆抗体

中文名称：兔抗 UFC1 多克隆抗体

英文名称：Anti-UFC1 rabbit polyclonal antibody

别名：HSPC155

抗原：UFC1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	UFC1 is an E2-like conjugating enzyme for ubiquitin-fold modifier-1 (UFM1; MIM 610553) (Komatsu et al., 2004 [PubMed 15071506]).
Applications:	ELISA, WB, IHC
Name of antibody:	UFC1
Immunogen:	Full length fusion protein
Full name:	ubiquitin-fold modifier conjugating enzyme 1
Synonyms:	HSPC155
SwissProt:	Q9Y3C8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and human cervical cancer

IHC Recommend dilution:	30-150
WB Predicted band size:	19 kDa
WB Positive control:	RAW264.7 cell lysate
WB Recommended dilution:	500-2000



