

兔抗 TDGF1 多克隆抗体

中文名称：兔抗 TDGF1 多克隆抗体

英文名称：Anti-TDGF1 rabbit polyclonal antibody

别名：CR; CRGF; CRIPTO

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TDGF1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with fore brain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	TDGF1
Immunogen:	Fusion protein of human TDGF1
Full name:	teratocarcinoma-derived growth factor 1
Synonyms:	CR; CRGF; CRIPTO

SwissProt:	P13385
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	40-200

