

## 兔抗 UHRF1BP1 多克隆抗体

中文名称: 兔抗 UHRF1BP1 多克隆抗体

英文名称: Anti-UHRF1BP1 rabbit polyclonal antibody

别名: ICBP90; C6orf107; dj349A12.1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: UHRF1BP1

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	UHRF1BP1 (UHRF1 binding protein 1) is a 1,440 amino acid protein that interacts with UHRF1 and may act as a negative regulator of cell growth. The UHRF1BP1 protein has been identified as a possible risk loci for systemic lupus erythematosus. The UHRF1BP1 gene is conserved in chimpanzee, dog, cow, mouse, rat, chicken, zebrafish and C.elegans, and maps to human chromosome 6p21.31.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UHRF1BP1
<b>Immunogen:</b>	Synthetic peptide of human UHRF1BP1
<b>Full name:</b>	UHRF1 binding protein 1

<b>Synonyms :</b>	ICBP90; C6orf107; dJ349A12.1
<b>SwissProt:</b>	Q6BDS2
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100

