

## 兔抗 UNC5B 多克隆抗体

中文名称：兔抗 UNC5B 多克隆抗体

英文名称：Anti-UNC5B rabbit polyclonal antibody

级 别：UNC5H2; p53RDL1

抗 原：UNC5B

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the netrin family of receptors. This particular protein mediates the repulsive effect of netrin-1 and is a vascular netrin receptor. This encoded protein is also in a group of proteins called dependence receptors (DpRs) which are involved in pro- and anti-apoptotic processes. Many DpRs are involved in embryogenesis and in cancer progression. Two alternatively spliced transcript variants have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UNC5B
<b>Immunogen:</b>	Synthetic peptide of human UNC5B

<b>Full name:</b>	unc-5 netrin receptor B
<b>Synonyms:</b>	UNC5H2; p53RDL1
<b>SwissProt:</b>	Q8IZJ1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100

