

兔抗 UQCRH 多克隆抗体

中文名称: 兔抗 UQCRH 多克隆抗体

英文名称: Anti-UQCRH rabbit polyclonal antibody

别名: QCR6; UQCR8

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: UQCRH

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondria I respiratory chain. This protein may mediate for mation of the complex between cytochromes c a nd c1.
Applications:	ELISA, WB, IHC
Name of antibody:	UQCRH
Immunogen:	Fusion protein of human UQCRH
Full name:	ubiquinol-cytochrome c reductase hinge protein
Synonyms:	QCR6; UQCR8

SwissProt:	P07919
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human prostate cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	11 kDa
WB Positive control:	Rat heart tissue and Mouse heart tissue lysates
WB Recommended dilution:	200-1000





