

兔抗 VPREB1 多克隆抗体

中文名称：兔抗 VPREB1 多克隆抗体

英文名称：Anti-VPREB1 rabbit polyclonal antibody

级别：IGI; IGVPB; VPREB; CD179a

抗原：VPREB1

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the immunoglobulin superfamily and is expressed selectively at the early stages of B cell development, namely, in pro B and early preB cells. This gene encodes the iota polypeptide chain that is associated with the Ig-mu chain to form a molecular complex which is expressed on the surface of pre-B cells. The complex is thought to regulate Ig gene rearrangements in the early steps of B-cell differentiation. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	VPREB1

Immunogen:	Synthetic peptide of human VPRED1
Full name:	pre-B lymphocyte 1
Synonyms:	IGI; IGVPB; VPRED; CD179a
SwissProt:	P12018
ELISA Recommended dilution:	500-1000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	20-100

