

兔抗 VPS36 多克隆抗体

中文名称: 兔抗 VPS36 多克隆抗体

英文名称: Anti-VPS36 rabbit polyclonal antibody

别名: EAP45; C13orf9; CGI-145

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: VPS36

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II (ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A similar protein complex in rat is associated with RNA polymerase elongation factor II.
Applications:	ELISA, IHC
Name of antibody:	VPS36
Immunogen:	Synthetic peptide of human VPS36
Full name:	vacuolar protein sorting 36 homolog
Synonyms:	EAP45; C13orf9; CGI-145

SwissProt:	Q86VN1
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human breast cancer and Human prostate cancer
IHC Recommend dilution:	30-150

