

## 兔抗 WSB1 多克隆抗体

中文名称:	兔抗 WSB1 多克隆抗体
英文名称:	Anti-WSB1 rabbit polyclonal antibody
别名:	SWIP1; WSB-1
储存:	冷冻 (-20℃)
抗原:	WSB1
宿主:	Rabbit
反应种属:	Human, Mouse
标记物:	Unconjugate
相关类别:	一抗
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the WD-protein sub family. This protein shares a high sequence identity to mouse and chick proteins. It contains several WD-repeats spanning most of the protein and an SOCS box in the C-terminus. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	WSB1
<b>Immunogen:</b>	Fusion protein of human WSB1
<b>Full name:</b>	WD repeat and SOCS box containing 1
<b>Synonyms :</b>	SWIP1; WSB-1

<b>SwissProt:</b>	Q9Y6I7
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	47 kDa
<b>WB Positive control:</b>	TM4 and NIH/3T3 cell
<b>WB Recommended dilution:</b>	500-2000

