

兔抗 ZBTB8A 多克隆抗体

中文名称: 兔抗 ZBTB8A 多克隆抗体

英文名称: Anti-ZBTB8A rabbit polyclonal antibody

别名: zinc finger and BTB domain containing 8A; BOZF1; ZBTB8; ZNF916A

储存: 冷冻 (-20℃)

抗原: ZBTB8A

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, WB, IHC
Name of antibody:	ZBTB8A
Immunogen:	Fusion protein of human ZBTB8A
Full name:	zinc finger and BTB domain containing 8A
Synonyms:	BOZF1; ZBTB8; ZNF916A
SwissProt:	Q96BR9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	50 kDa
WB Positive control:	293T cell lysate

WB Recommended dilution:

500-2000



