

兔抗 TEX11 多克隆抗体

中文名称: 兔抗 TEX11 多克隆抗体

英文名称: Anti-TEX11 rabbit polyclonal antibody

别名: testis expressed 11; TGC1; TSGA3; SPGFX2

储存: 冷冻 (-20℃)

抗原: TEX11

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TEX11
Immunogen:	Fusion protein of human TEX11
Full name:	testis expressed 11
Synonyms:	TGC1; TSGA3; SPGFX2
SwissProt:	Q8IYF3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human esophagus cancer

IHC Recommend dilution:

50-300



