

兔抗 TGM4 多克隆抗体

中文名称：兔抗 TGM4 多克隆抗体

英文名称：Anti-TGM4 rabbit polyclonal antibody

别名：TGP; hTGP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TGM4

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Associated with the mammalian reproductive process. Catalyzes the cross-linking of proteins and the conjugation of polyamines to specific proteins in the seminal tract.
Applications:	ELISA, IHC
Name of antibody:	TGM4
Immunogen:	Synthetic peptide of human TGM4
Full name:	transglutaminase 4
Synonyms :	TGP; hTGP
SwissProt:	P49221
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human prostate cancer

IHC Recommend dilution:

10-50

