

兔抗 THAP12 多克隆抗体

- 中文名称: 兔抗 THAP12 多克隆抗体
- 英文名称: Anti-THAP12 rabbit polyclonal antibody
- 别名: THAP domain containing 12; DAP4; THAP0; PRKRIR; P52rIPK
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 抗原: THAP12
- 宿主: Rabbit
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.
Applications:	ELISA, IHC
Name of antibody:	THAP12
Immunogen:	Fusion protein of human THAP12
Full name:	THAP domain containing 12
Synonyms:	DAP4; THAP0; PRKRIR; P52rIPK
SwissProt:	O43422
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer

IHC Recommend dilution:

30-150

