

兔抗 PARD3 多克隆抗体

中文名称：兔抗 PARD3 多克隆抗体

英文名称： Anti-PARD3 rabbit polyclonal antibody

别名： Baz; ASIP; PAR3; PARD-3; PARD3A; SE2-5T2; PPP1R118; SE2-5L16; SE2-5LT1; PAR3alpha

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PARD3

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene.
Applications:	WB, IF
Name of antibody:	PARD3
Immunogen:	Synthesized peptide derived from internal of human PARD3 .
Full name:	par-3 family cell polarity regulator

Synonyms :	Baz; ASIP; PAR3; PARD-3; PARD3A; SE2-5T2; PPP1R118; SE2-5L16; SE2-5LT1; PAR3alpha
SwissProt:	Q8TEW0
WB Predicted band size:	151 kDa
WB Positive control:	COLO205 cells lysate
WB Recommended dilution:	500-3000
IF positive control:	HepG2 cells
IF Recommend dilution:	100-500



