

## 兔抗 THRA 多克隆抗体

中文名称：兔抗 THRA 多克隆抗体

英文名称：Anti-THRA rabbit polyclonal antibody

别名：AR7; EAR7; ERBA; CHNG6; ERBA1; NR1A1; THRA1; THRA2; ERB-T-1; c-ERBA-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：THRA/THRB

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a nuclear hormone receptor for triiodothyronine. It is one of the several receptors for thyroid hormone, and has been shown to mediate the biological activities of thyroid hormone. Knockout studies in mice suggest that the different receptors, while having certain extent of redundancy, may mediate different functions of thyroid hormone. Alternatively spliced transcript variants encoding distinct isoforms have been reported.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	THRA/THRB
<b>Immunogen:</b>	Synthetic peptide of human THRA/THRB

<b>Full name:</b>	thyroid hormone receptor, alpha/beta
<b>Synonyms:</b>	AR7; EAR7; ERBA; CHNG6; ERBA1; NR1A1; THRA1; THRA2; ERB-T-1; c-ERBA-1; GRTH; PRTH; THR1; ERBA2; NR1A2; THRB1; THRB2; C-ERBA-2; C-ERBA-BETA
<b>SwissProt:</b>	P10827/P10828
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	55 kDa
<b>WB Positive control:</b>	Human fetal brain tissue lysate
<b>WB Recommended dilution:</b>	500-2000

