

兔抗 TLE1 多克隆抗体

中文名称：兔抗 TLE1 多克隆抗体

英文名称：Anti-TLE1 rabbit polyclonal antibody

别名：transducin like enhancer of split 1; ESG; ESG1; GRG1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TLE1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Transcriptional corepressor that binds to a number of transcription factors. Inhibits NF-kappa-B-regulated gene expression. Inhibits the transcriptional activation mediated by FOXA2, and by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Unusual function as coactivator for ESRRG.
Applications:	ELISA, WB, IHC
Name of antibody:	TLE1
Immunogen:	Fusion protein of human TLE1
Full name:	transducin like enhancer of split 1

Synonyms:	ESG; ESG1; GRG1
SwissProt:	Q04724
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	83 kDa
WB Positive control:	Human fetal liver tissue lysate
WB Recommended dilution:	200-1000



