

## 兔抗 PCID2 多克隆抗体

中文名称：兔抗 PCID2 多克隆抗体

英文名称：Anti-PCID2 rabbit polyclonal antibody

别名：F10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

抗原：PCID2

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a component of the TREX-2 complex (transcription and export complex 2), which regulates mRNA export from the nucleus. This protein regulates expression of Mad2 mitotic arrest deficient-like 1, a cell division checkpoint protein. This protein also interacts with and stabilizes Brca2 (breast cancer 2) protein. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PCID2
<b>Immunogen:</b>	Fusion protein of human PCID2
<b>Full name:</b>	PCI domain containing 2
<b>Synonyms:</b>	F10
<b>SwissProt:</b>	Q5JVF3

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Mouse testis tissue, Rat heart tissue lysates
<b>WB Recommended dilution:</b>	500-2000



