

兔抗 TMA16 多克隆抗体

中文名称: 兔抗 TMA16 多克隆抗体

英文名称: Anti-TMA16 rabbit polyclonal antibody

别名: translation machinery associated 16 homolog; C4orf43

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TMA16

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	TMA16
Immunogen:	Fusion protein of human TMA16
Full name:	translation machinery associated 16 homolog
Synonyms:	C4orf43
SwissProt:	Q96EY4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human gastric cancer
IHC Recommend dilution:	100-200

