

兔抗 PDCD2 多克隆抗体

中文名称：兔抗 PDCD2 多克隆抗体

英文名称：Anti-PDCD2 rabbit polyclonal antibody

别名：RP8, ZMYND7

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PDCD2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear protein expressed in a variety of tissues. Expression of this gene has been shown to be repressed by B-cell CLL/lymphoma 6 (BCL6), a transcriptional repressor required for lymph node germinal center development, suggesting that BCL6 regulates apoptosis by its effects on this protein. Alternative splicing results in multiple transcript variants and pseudogenes have been identified on chromosomes 9 and 12.
Applications:	ELISA, IHC
Name of antibody:	PDCD2
Immunogen:	Synthetic peptide of human PDCD2
Full name:	Programmed cell death 2

Synonyms :	RP8, ZMYND7
SwissProt:	Q16342
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human ovarian cancer and Human tonsil
IHC Recommend dilution:	50-200



