

## 兔抗 PDCD4 多克隆抗体

中文名称: 兔抗 PDCD4 多克隆抗体

英文名称: Anti-PDCD4 rabbit polyclonal antibody

别名: programmed cell death 4; H731

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PDCD4

反应种属: Human, Mouse, Rat



标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding . Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PDCD4
<b>Immunogen:</b>	Fusion protein of human PDCD4
<b>Full name:</b>	programmed cell death 4
<b>Synonyms:</b>	H731
<b>SwissProt:</b>	Q53EL6
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human ovarian cancer and Human tonsil
<b>IHC Recommend dilution:</b>	40-200



**mlbio** 梅联生物  
Good elisakit producers