

兔抗 PDCD7 多克隆抗体

中文名称: 兔抗 PDCD7 多克隆抗体

英文名称: Anti-PDCD7 rabbit polyclonal antibody

别名: ES18, HES18

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PDCD7

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a 59 kDa protein that is associated with the U11 small nuclear ribonucleoprotein (snRNP), which is a component of the minor U12-type spliceosome responsible for catalyzing pre-mRNA splicing of U12-type introns.
Applications:	ELISA, IHC
Name of antibody:	PDCD7
Immunogen:	Synthetic peptide of human PDCD7
Full name:	Programmed cell death 7
Synonyms :	ES18, HES18
SwissProt:	Q8N8D1
ELISA Recommended dilution:	2000-10000

IHC positive control:	Human colon cancer and Human lung cancer
IHC Recommend dilution:	50-200

