

兔抗 PDRG1 多克隆抗体

中文名称: 兔抗 PDRG1 多克隆抗体

英文名称: Anti-PDRG1 rabbit polyclonal antibody

别名: p53 and DNA damage regulated 1; PDRG; C20orf126

抗原: PDRG1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:	May play a role in chaperone-mediated protein folding.
Applications:	ELISA, IHC
Name of antibody:	PDRG1
Immunogen:	Fusion protein of human PDRG1
Full name:	p53 and DNA damage regulated 1
Synonyms:	PDRG; C20orf126
SwissProt:	Q9NUG6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

