

兔抗 PDXDC1 多克隆抗体

中文名称: 兔抗 PDXDC1 多克隆抗体

英文名称: Anti-PDXDC1 rabbit polyclonal antibody

别名: LP8165

储存: 冷冻 (-20℃) 避光

抗原: PDXDC1

宿主: Rabbit

反应种属: Human Mouse Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Belongs to the group II decarboxylase family. 2 isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	PDXDC1
Immunogen:	Synthesized peptide derived from internal of human PDXDC1.
Full name:	pyridoxal-dependent decarboxylase domain containing 1
Synonyms :	LP8165
SwissProt:	Q6P996
WB Predicted band size:	87 kDa
WB Positive control:	HepG2 cells and Huvec cells lysates
WB Recommended dilution:	500-3000

