

兔抗 PER3 多克隆抗体

中文名称：兔抗 PER3 多克隆抗体

英文名称：Anti-PER3 rabbit polyclonal antibody

别名：period circadian regulator 3; GIG13; FASPS3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PER3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been linked to sleep disorders. Multiple transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, IHC
Name of antibody:	PER3
Immunogen:	Fusion protein of human PER3
Full name:	period circadian regulator 3
Synonyms:	GIG13; FASPS3
SwissProt:	P56645
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human cervical cancer
IHC Recommend dilution:	50-300



