

兔抗 PGR 多克隆抗体

中文名称：兔抗 PGR 多克隆抗体

英文名称：Anti-PGR rabbit polyclonal antibody

别名：PR, NR3C3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PGR

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the steroid receptor superfamily. The encoded protein mediates the physiological effects of progesterone, which plays a central role in reproductive events associated with the establishment and maintenance of pregnancy. This gene uses two distinct promoters and translation start sites in the first exon to produce two isoforms, A and B. The two isoforms are identical except for the additional 165 amino acids found in the N-terminus of isoform B and mediate their own response genes and physiologic effects with little overlap.
Applications:	ELISA, IHC
Name of antibody:	PGR

Immunogen:	Synthetic peptide of human PGR
Full name:	Progesterone receptor
Synonyms :	PR, NR3C3
SwissProt:	P06401
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human breast cancer and Human cervical cancer
IHC Recommend dilution:	25-100



