

兔抗 PISD 多克隆抗体

中文名称: 兔抗 PISD 多克隆抗体

英文名称: Anti-PISD rabbit polyclonal antibody

别名: phosphatidylserine decarboxylase; PSD; PSDC; PSSC; DJ858B16; dj858B16.2

储存: 冷冻 (-20℃)

抗原: PISD

宿主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene catalyzes the conversion of phosphatidylserine to phosphatidylethanolamine in the inner mitochondrial membrane. The encoded protein is active in phospholipid metabolism and inter organelle trafficking of phosphatidylserine.
Applications:	ELISA, IHC
Name of antibody:	PISD
Immunogen:	Fusion protein of human PISD
Full name:	phosphatidylserine decarboxylase
Synonyms:	PSD; PSDC; PSSC; DJ858B16; dj858B16.2
SwissProt:	Q9UG56
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300

