

兔抗 PLA2G15 多克隆抗体

中文名称：兔抗 PLA2G15 多克隆抗体

英文名称： Anti-PLA2G15 rabbit polyclonal antibody

别名： phospholipase A2 group XV; ACS; LLPL; LPLA2; LYPLA3; GXVPLA2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PLA2G15

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. The protein encoded by this gene hydrolyzes lysophosphatidylcholine to glycerophosphorylcholine and a free fatty acid. This enzyme is present in the plasma and thought to be associated with high-density lipoprotein. A later paper contradicts the function of this gene. It demonstrates that this gene encodes a lysosomal enzyme instead of a lysophospholipase and has both calcium-independent phospholipase A2 and transacylase activities.
Applications:	ELISA, IHC
Name of antibody:	PLA2G15

Immunogen:	Fusion protein of human PLA2G15
Full name:	phospholipase A2 group XV
Synonyms:	ACS; LLPL; LPLA2; LYPLA3; GXVPLA2
SwissProt:	Q8NCC3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human prostate cancer
IHC Recommend dilution:	50-200



