

兔抗 PLK4 多克隆抗体

中文名称：兔抗 PLK4 多克隆抗体

英文名称：Anti-PLK4 rabbit polyclonal antibody

别名：SAK; STK18; MCCR2

抗原：PLK4

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the polo family of serine/threonine protein kinases. The protein localizes to centrioles, complex microtubule-based structures found in centrosomes, and regulates centriole duplication during the cell cycle. Three alternatively spliced transcript variants that encode different protein isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	PLK4
Immunogen:	Fusion protein of human PLK4
Full name:	polo-like kinase 4
Synonyms:	SAK; STK18; MCCR2

SwissProt:	O00444
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	35-200

