

兔抗 PMP22 多克隆抗体

中文名称：兔抗 PMP22 多克隆抗体

英文名称：Anti-PMP22 rabbit polyclonal antibody

别名：peripheral myelin protein 22; DSS; GAS3; HNPP; CMT1A; CMT1E; GAS-3; Sp110; HMSNIA

相关类别：一抗

抗原：PMP22

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an integral membrane protein that is a major component of myelin in the peripheral nervous system. Studies suggest two alternately used promoters drive tissue-specific expression. Various mutations of this gene are causes of Charcot-Marie-Tooth disease Type IA, Dejerine-Sottas syndrome, and hereditary neuropathy with liability to pressure palsies. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	PMP22

Immunogen:	Synthetic peptide of human PMP22
Full name:	peripheral myelin protein 22
Synonyms:	DSS; GAS3; HNPP; CMT1A; CMT1E; GAS-3; Sp110; HMSNIA
SwissProt:	Q01453
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	18 kDa
WB Positive control:	Human cerebrum tissue , A172 cell , Mouse lung tissue , 293 cell , 231 cell , Mouse liver tissue , Mouse brain tissue lysate
WB Recommended dilution:	200-1000



