

兔抗 POLA1 多克隆抗体

中文名称: 兔抗 POLA1 多克隆抗体

英文名称: Anti-POLA1 rabbit polyclonal antibody

别名: NSX; POLA; p180

储存: 冷冻 (-20℃)

抗原: POLA1

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the catalytic subunit of DNA polymerase α , which together with a regulatory and two primase subunits, forms the DNA polymerase alpha complex. The catalytic subunit plays an essential role in the initiation of DNA replication.
Applications:	ELISA, IHC
Name of antibody:	POLA1
Immunogen:	Synthetic peptide of human POLA1
Full name:	polymerase (DNA) alpha 1, catalytic subunit
Synonyms:	NSX; POLA; p180

SwissProt:	P09884
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	20-100

