

## 兔抗 TP53 多克隆抗体

中文名称：兔抗 TP53 多克隆抗体

英文名称：Anti-TP53 rabbit polyclonal antibody

别名：P53, BCC7, LFS1, TRP53

相关类别：一抗

抗原：TP53

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, there by inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons from identi

	cal transcript variants.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	TP53
<b>Immunogen:</b>	Fusion protein of human TP53
<b>Full name:</b>	tumor protein p53
<b>Synonyms:</b>	P53; BCC7; LFS1; TRP53
<b>SwissProt:</b>	P04637
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	44 kDa
<b>WB Positive control:</b>	293T cell lysate
<b>WB Recommended dilution:</b>	500-2000

