

兔抗 PPIL1 多克隆抗体

中文名称: 兔抗 PPIL1 多克隆抗体

英文名称: Anti-PPIL1 rabbit polyclonal antibody

别名: CYPL1; hCyPX; PPIase; CGI-124

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PPIL1

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene is a member of the cyclophilin family of peptidylprolyl isomerases (PPIases). The cyclophilins are a highly conserved, ubiquitous family, members of which play an important role in protein folding, immunosuppression by cyclosporin A, and infection of HIV-1 virions. Based on similarity to other PPIases, this protein could accelerate the folding of proteins and might catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.
Applications:	ELISA, WB, IHC
Name of antibody:	PPIL1
Immunogen:	Full length fusion protein

Full name:	peptidylprolyl isomerase (cyclophilin)-like 1
Synonyms :	CYPL1; hCyPX; PPIase; CGI-124
SwissProt:	Q9Y3C6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human prostate cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	18 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	500-2000

