

兔抗 TRARG1 多克隆抗体

中文名称: 兔抗 TRARG1 多克隆抗体

英文名称: Anti-TRARG1 rabbit polyclonal antibody

别名: trafficking regulator of GLUT4 (SLC2A4) 1; BEC-1;

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TRARG1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Regulates insulin-mediated adipose tissue glucose uptake and transport by modulation of SLC2A4 recycling. Not required for SLC2A4 membrane fusion upon an initial stimulus, but rather is necessary for proper protein recycling during prolonged insulin stimulation.
Applications:	ELISA, IHC
Name of antibody:	TRARG1
Immunogen:	Synthetic peptide of human TRARG1
Full name:	trafficking regulator of GLUT4 (SLC2A4) 1
Synonyms:	BEC-1; DSPB1; LOST1; TUSC5; IFITMD3
SwissProt:	Q8IXB3
IHC positive control:	Human thyroid cancer

IHC Recommend dilution:	30-150
-------------------------	--------

