

兔抗 PSD4 多克隆抗体

中文名称:	兔抗 PSD4 多克隆抗体
英文名称:	Anti-PSD4 rabbit polyclonal antibody
别名:	pleckstrin and Sec7 domain containing 4; TIC; EFA6B
储存:	冷冻 (-20℃)
抗原:	PSD4
宿主:	Rabbit
反应种属:	Human, Mouse
相关类别:	一抗
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	Guanine nucleotide exchange factor for ARF6 and ARL14 /ARF7. Through ARL14 activation, controls the movement of MHC class II-containing vesicles along the actin cytoskeleton in dendritic cells. Involved in membrane recycling. Interacts with several phosphatidylinositol phosphate species, including phosphatidylinositol 3,4-bisphosphate, phosphatidylinositol 3,5-bisphosphate and phosphatidylinositol 4,5-bisphosphate.
Applications:	ELISA, IHC
Name of antibody:	PSD4
Immunogen:	Fusion protein of human PSD4

Full name:	pleckstrin and Sec7 domain containing 4
Synonyms:	TIC; EFA6B
SwissProt:	Q8NDX1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150

