

兔抗 TRIM63 多克隆抗体

中文名称: 兔抗 TRIM63 多克隆抗体

英文名称: Anti-TRIM63 rabbit polyclonal antibody

别名: IRF, SMRZ, MURF1, MURF2, RNF28

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TRIM63

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Synthetic peptide of human TRIM63
Applications:	ELISA, WB, IHC
Name of antibody:	TRIM63
Immunogen:	Synthetic peptide of human TRIM63
Full name:	tripartite motif containing 63, E3 ubiquitin protein ligase
Synonyms:	IRF, SMRZ, MURF1, MURF2, RNF28
SwissProt:	Q969Q1
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer;

IHC Recommend dilution:	50-200
WB Predicted band size:	40 kDa
WB Positive control:	NIH/3T3 cells and human fetal muscle tissue
WB Recommended dilution:	1000-5000



