

兔抗 SEPTIN1 多克隆抗体

中文名称：兔抗 SEPTIN1 多克隆抗体

英文名称：Anti-SEPTIN1 rabbit polyclonal antibody

别名：septin 1; LARP; SEP1; DIFF6; SEPT1; PNUTL3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SEPTIN1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis and the maintenance of cellular morphology. This gene encodes a protein that can form homo- and heterooligomeric filaments, and may contribute to the formation of neurofibrillary tangles in Alzheimer's disease. Alternatively spliced transcript variants have been found but the full-length nature of these variants has not been determined. [provided by RefSeq, Dec 2012]
Applications:	ELISA, IHC
Name of antibody:	SEPTIN1
Immunogen:	Fusion protein of human SEPTIN1
Full name:	septin 1
Synonyms:	LARP; SEP1; DIFF6; SEPT1; PNUTL3

SwissProt:	Q8WYJ6
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	40-200

