

## 兔抗 SHB 多克隆抗体

中文名称：兔抗 SHB 多克隆抗体

英文名称： Anti-SHB rabbit polyclonal antibody

别名： bA3J10.2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SHB

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	Adapter protein which regulates several signal transduction cascades by linking activated receptors to downstream signaling components. May play a role in angiogenesis by regulating FGFR1, VEGFR2 and PDGFR signaling. May also play a role in T-cell antigen receptor/TCR signaling, interleukin-2 signaling, apoptosis and neuronal cells differentiation by mediating basic-FGF and NGF-induced signaling cascades. May also regulate IRS1 and IRS2 signaling in insulin-producing cells.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	SHB

<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human SHB.
<b>Full name:</b>	Src homology 2 domain containing adaptor protein B
<b>Synonyms :</b>	bA3J10.2
<b>SwissProt:</b>	Q15464
<b>WB Predicted band size:</b>	55 kDa
<b>WB Positive control:</b>	Huvec cells lysate
<b>WB Recommended dilution:</b>	500-3000

