

兔抗 ADAM2 多克隆抗体

中文名称：兔抗 ADAM2 多克隆抗体

英文名称：Anti-ADAM2 rabbit polyclonal antibody

别名：CT15; FTNB; PH30; CRYN1; CRYN2; PH-30b

相关类别：一抗

储存：冷冻（-20℃）

抗原：ADAM2

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. This member is a subunit of an integral sperm membrane glycoprotein called fertilin, which plays an important role in sperm-egg interactions.
Applications:	ELISA, IHC
Name of antibody:	ADAM2

Immunogen:	Fusion protein of human ADAM2
Full name:	ADAM metallopeptidase domain 2
Synonyms :	CT15; FTNB; PH30; CRYN1; CRYN2; PH-30b
SwissProt:	Q6P2G0
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	15-50

