

兔抗 ADAP2 多克隆抗体

中文名称：兔抗 ADAP2 多克隆抗体

英文名称：Anti-ADAP2 rabbit polyclonal antibody

别名：ArfGAP with dual PH domains 2; CENTA2; cent-b; HSA272195

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ADAP2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	GTPase-activating protein for the ADP ribosylat ion factor family (Potential). Binds phosphatidyl inositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4). Possesses a stoichiometry of two binding sites for InsP4 with identical affinity.
Applications:	ELISA, IHC
Name of antibody:	ADAP2
Immunogen:	Fusion protein of human ADAP2
Full name:	ArfGAP with dual PH domains 2
Synonyms:	CENTA2; cent-b; HSA272195
SwissProt:	Q9NPF8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

