

兔抗 ADGRF1 多克隆抗体

中文名称: 兔抗 ADGRF1 多克隆抗体

英文名称: Anti-ADGRF1 rabbit polyclonal antibody

别名: adhesion G protein-coupled receptor F1; PGR19; GPR110; KPG_012; hGPCR36

抗原: ADGRF1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Orphan receptor.
Applications:	ELISA, IHC
Name of antibody:	ADGRF1
Immunogen:	Synthetic peptide of human ADGRF1
Full name:	adhesion G protein-coupled receptor F1
Synonyms:	PGR19; GPR110; KPG_012; hGPCR36
SwissProt:	Q5T601
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	40-200

