

## 兔抗 ADGRF5 多克隆抗体

中文名称: 兔抗 ADGRF5 多克隆抗体

英文名称: Anti-ADGRF5 rabbit polyclonal antibody

别名: adhesion G protein-coupled receptor F5; GPR116; KPG\_001

储存: 冷冻 (-20℃)

抗原: ADGRF5

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

<b>Background:</b>	Receptor that plays a critical role in lung surfactant homeostasis. May play a role in controlling adipocyte function.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ADGRF5
<b>Immunogen:</b>	Synthetic peptide of human ADGRF5
<b>Full name:</b>	adhesion G protein-coupled receptor F5
<b>Synonyms:</b>	GPR116; KPG_001
<b>SwissProt:</b>	Q8IZF2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	40-200

