

兔抗 ADORA2B 多克隆抗体

中文名称：兔抗 ADORA2B 多克隆抗体

英文名称：Anti-ADORA2B rabbit polyclonal antibody

级别：ADORA2

抗原：ADORA2B

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an adenosine receptor that is a member of the G protein-coupled receptor superfamily. This integral membrane protein stimulates adenylate cyclase activity in the presence of adenosine. This protein also interacts with netrin-1, which is involved in axon elongation. The gene is located near the Smith-Magenis syndrome region on chromosome 17.
Applications:	ELISA, IHC
Name of antibody:	ADORA2B
Immunogen:	Synthetic peptide of human ADORA2B
Full name:	adenosine A2b receptor
Synonyms:	ADORA2

SwissProt:	P29275
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

