

## 兔抗 AHCYL1 多克隆抗体

中文名称：兔抗 AHCYL1 多克隆抗体

英文名称：Anti-AHCYL1 rabbit polyclonal antibody

别名：adenosylhomocysteinase like 1; DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AHCYL1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene interacts with inositol 1,4,5-trisphosphate receptor, type 1 and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding two different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	AHCYL1
<b>Immunogen:</b>	Fusion protein of human AHCYL1
<b>Full name:</b>	adenosylhomocysteinase like 1
<b>Synonyms:</b>	DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA

<b>SwissProt:</b>	O43865
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	50-300

